

POWERSTROLL Range

LIGHTWEIGHT POWERSTROLL



This excellent lightweight Powerstroll will turn most manual wheelchairs into carer controlled powered wheelchairs.

Product features

- Allows speed of up to 4.5kmph*
- Travels up to 16km (10 miles) on a full battery charge*
- Fits most wheelchairs with seat widths from 14" to 20" (including narrow track wheelchairs such as the Lightweight Aluminium)
- Fits 2.2cm-2.5cm (7/8"-1") diameter tubing
- Uses existing wheelchair handgrips
- Lightweight at 14kg and less than 8kg (17lb) with battery pack removed for easy transportation
- Trigger handset with variable speed dial and direction button to suit the attendant and user requirements
- The Powerstroll can be easily fitted and removed
- Twin 'Drive' wheels with 8" black non-marking PU tyres
- Removable fabric battery pack for easy storage and transportation
- Uses 12v 20ah battery (airline friendly)
- Lightweight 100w 12 DC motor making the unit very portable
- Easy to use lift up strap to increase manoeuvrability.
- Weight capacity 115kg (18st)
- Adaptor brackets are available to allow the telescopic bar at the front of the powerpack to be installed at its standard height or at a new height approx. 2" (5cm) lower
- This is ideal for installation on wheelchairs, that have a higher or lower horizontal bar or lack a hemi-shape (such as the XS2 and SD2)
- The brackets are supplied in pairs. No further fasteners are required

CODE: PWCPP010AU



| Code | Weight Capacity | Total Weight | Weight Without Battery | Weight Battery Pack Only | Maximum Gradient | Maximum Range |
|------------|-----------------|--------------|------------------------|--------------------------|------------------|----------------|
| PWCPP010AU | 115kg (18st) | 14kg (31lb) | 8kg (17lb) | 7kg (14lb) | 6° | 16km (10miles) |

www.drive-medical.com.au

drive

DeVilbiss
HEALTHCARE

POWERSTROLL Range



S DRIVE POWERSTROLL

Superb Comfort, Style & Performance...

- Allows speed of up to 4mph*
- Travels up to 10 miles (16 km) on a full battery charge*
- The Powerstroll can easily be fitted or removed within seconds
- The excellent P G Drives S-Drive controller with fully programmability & comprehensive fault management system
- New push and slide battery removal system
- It is very portable, making transportation simple
- Fits most standard wheelchairs up to 20" seat width
- Hand control with variable speed dial, forward/reverse function, illuminated battery gauge and fingertip trigger control
- Unique one-click docking station for ease of assembly (patented)
- Comes complete with a Carry Holdall with shoulder strap
- Simple to manoeuvre and reverse
- The Powerstroll has a powerful 150 watt (200 watt-HD), 24V motor, capable of coping with most terrains
- Maintenance free airline friendly 24V 12ah valve regulated sealed lead acid battery pack
- 2A charger as standard for off-board charging
- Weight capacity 135kg (21st), this includes the wheelchair and user combined
- Heavy Duty model - weight capacity 180kg (28st)

It takes the strain out of pushing a heavy wheelchair and makes easy work of hills, ramps and uneven surfaces. Its unique one click docking station (patented) means that it can be fitted or removed easily in seconds.

* Subject to conditions of use

CODE: PWCPP009AU
PWCPP009HDA – Heavy Duty



| Code | Weight Capacity | Total Weight | Weight Without Battery | Weight Battery Pack Only | Maximum Gradient | Maximum Range |
|--------------|-----------------|---------------|------------------------|--------------------------|------------------|----------------|
| PWCPP009AU | 135kg (21st) | 19.2kg (42lb) | 10.2kg (22lb) | 9kg (20lb) | 8° | 16km (10miles) |
| PWCPP009HDAU | 180kg (28st) | 19.2kg (42lb) | 10.2kg (22lb) | 9kg (20lb) | 8° | 16km (10miles) |